DUR-LAB-Sil®

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Page 3 - en: Instructions for use

Page 4 - fr: Manuel d'instruction

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Strana 7 - bs: Uputstvo za upotrebu

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Instructions for use

Patient target group

None, as DUR-LAB-Sil does not come into contact with the human body.

Intended user

Dental technician

Composition

Qualitative data:

Base paste:

Polydimethylsiloxanes with terminal hydroxyl groups, quartz powder, silicon dioxide

Catalyst paste:

Tetraethoxysilane, tin catalyst, silicone oil

Quantitative data in % w/w:

Base paste:

Polydimethylsiloxanes with terminal hydroxyl groups 30, quartz powder 60, silicon dioxide 10

Catalyst paste:

Tetraethoxysilane 49, tin catalyst 1, silicone oil 50

Intended purpose of the product

DUR-LAB-Sil is a kneadable silicone based on condensation-curing silicones for laboratory applications such as bite walls and blockouts.

Indications

For making bite walls and blocking out undercuts.

Contraindications

Do not use if allergic to the ingredients.

Performance characteristics of the product

DUR-LAB-Sil cures by mixing base paste with catalyst paste.

Specific precautions

None

Application

Processing with DUR-LAB-Sil paste:

Depending on the desired quantity, 1 strand of catalyst paste is dosed per 1 spoon of base mass per diameter of this dosing spoon. Knead both masses until the colour is uniform throughout. The processing time when used with DUR-LAB-Sil hardener paste is approx. 2–5 min. at room temperature. After 8 minutes the material should be set.

Working instructions

The working instructions and precautionary measures described in the package leaflet accompanying the product must be observed in every instance. If, despite having followed these instructions, undesirable events occur, we request that you inform us in every instance, describing the accompanying circumstances as precisely as possible.

Duration of application

For caries prophylaxis, treatment should be carried out 2–3 times a year (or 3–4 times if there is an increased susceptibility to caries).

Warnings**Residual risks**

H226	Flammable liquid and vapour.
H330	Can be fatal if inhaled.
H315	Causes skin irritation.
H318	Causes serious eye damage.
H317	May provoke allergic skin reactions.
H341	May cause genetic defects.
H360FD	May impair fertility. May harm the child in the womb
H371	May cause organ damage.
H373	May cause organ damage through prolonged or multiple exposure.
H411	Toxic to aquatic organisms, with long-lasting effect.

Danger Danger

Attention Warning

Side effects

Adverse side effects are not known.

Interactions

Interaction of the product and its metabolites with other medical devices, medicinal products or other substances is not known.

Storage

Minimum shelf life at storage under 23 °C is stated on label. Store sealed and under 25 °C!

Minimum shelf life

Do not use after the expiry date.

Variants

[REF] WP5060 DUR-LAB-Sil, 2600g

[REF] WP5300 DUR-LAB-Sil, 5000g

[REF] WP5067 DUR-LAB-Sil, Activator Paste, 40.0g

Disposal

Unused material should be disposed of appropriately after curing.

All serious incidents related to the product shall be reported to the manufacturer and to the competent authority of the Member State where the user and/or the patient is established.

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Willmann & Pein GmbH