WP Metalprime®

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Page 3 - en: Instructions for use

Page 4 - fr: Manuel d'instruction

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Strana 7 - bs: Uputstvo za upotrebu

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Instructions for use

Patient target group

Children, adolescents, adults and older patients

Intended users

Dental technician and dentist

Composition

Qualitative data:

Methacrylate, water, ethanol, MDP

Quantitative data in % G/G:

Water/Ethanol 97, MDP 2

Intended purpose of the product

Adhesion fluid for the bonding of metals and plastics/composites

Indications

WP Metalprime is used to improve the adhesion of plastics and composites to metals of all kinds.

Contraindications

WP Metalprime may be used only by dentists and dental technicians. It contains methacrylate. The product therefore should not be used on / by persons known to be allergic to methacrylate. If allergic reactions occur, contact/visit a doctor immediately.

Performance characteristics of the product

Adhesive bond of metal with plastic (methacrylate) by polymerisation.

Specific precautions

Avoid sparks and direct sunlight. Use product in well-ventilated workrooms.

Application

Sandblast the metal surfaces with aluminium oxide and then clean them in an air jet (recommended grain size 110 at 2.8 bar). Avoid contamination with release agents, oil or moisture.

Put 2-3 drops of the WP Metalprime into a suitable container. Apply a thin layer to the adhesive surface with the enclosed brush and allow to dry. If necessary, a second coat can also be applied. Apply the plastic to the metal immediately after the surface treatment. Do not use the enclosed brush for different products, but only for WP Metalprime.

Areas of application

- Repairs of denture clasps, metal frameworks and attachments.
Roughen the metal bonding surface with a carborundum sander, prepare the bonding surface as described above, and build up the autopolymer in the usual way.
- Opaque to metal
Composite crowns and composite veneers.
Prepare the surfaces as described above and apply opaque as usual. Polymerise the material according to the manufacturer's recommendations.
- Plastic cement to metal
Metal inlays, crowns and bridges.
Prepare the metal surface as described above and bond the restoration with resin cement.
- Connections made of plastics with orthodontic apparatuses made of metal
Prepare the bonding surface as described above and apply the orthodontic resin to the apparatus as usual.
- Denture base Magnet attachments
Prepare the surface as described above and attach as usual with autopolymer.
- Denture base
Plastic on metal frameworks, clasps or attachments.
Sandblast the bonding surface with aluminium oxide and clean it in an air jet. Make a wax model of the denture. Clean the bonding surface with superheated steam after the wax-out and remove the wax residues completely with a solvent. Then apply WP Metalprime and plug and polymerise the denture base in the usual way.

- Model casting/Veneer crowns

Blast retentions with aluminium oxide. Apply WP Metalprime with a clean brush and dry with a hot air gun or compressed air compressor (oil-free) with water separator.

Warnings**Residual risks**

The processing instructions and precautionary measures described in the instructions for use accompanying the product must always be observed.

If, despite following these instructions, an adverse reaction occurs in the patient or other adverse events occur, we request that you inform us in any case, describing the accompanying circumstances as precisely as possible.

Side effects

Adverse side effects are not known. Risks in case of overdose are not known.

Interactions

Interaction of the product and its metabolites with other medical devices, medicinal products or other substances is not known.

Storage

It is highly flammable and must be stored in a cool and dry place. WP Metalprime is volatile, therefore the container should be closed immediately after use

Minimum shelf life

Do not use after the expiry date.

Variants

REFWP5593 2 x 3,0 ml liquid + disposable brush and brush holder

Disposal

Unused material must be disposed of according to guidelines for methacrylate waste after the solvents (water/ethanol) have evaporated.

All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient is established.

We know how

dent a pharm
Produktionsgesellschaft mbH
Schusterring 35
25355 Barmstedt/Germany
Fon: +(49) 4123/9225-0
Fax: +(49) 4123/9225-48
www.dent-a-pharm.de
info@dent-a-pharm.de

Willmann & Pein GmbH
Schusterring 35
25355 Barmstedt/Germany
Fon: +(49) 4123/9228-0
Fax: +(49) 4123/9228-49
www.wp-dental.de
info@wp-dental.de

Made in Germany

H225 Highly flammable liquid and vapour
H314 Causes severe skin burns and eye damage

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