# LiWa®

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Page 3 - en: Instructions for use

Page 4 - fr: Manuel d'instruction

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Strana 7 - bs: Uputstvo za upotrebu

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## Instructions for use

### Patient target group

None, as LiWa does not come into contact with the human body.

### Intended user

Dental technician

### Composition

Qualitative data: Acrylates, photoinitiators, accelerators  
Quantitative data in % w/w: > 90 % acrylates

### Intended purpose of the product

LiWa is a light-curing modelling material with which metal parts (crowns, bridges) are made in the laboratory.

### Indications

For the production of plastic moulded parts for casting technology.

### Contraindications

When modelling, vapours may form due to heating. Do not use if allergic to the ingredients.

### Performance characteristics of the product

LiWa is a ready-to-use, light-curing modelling clay with a firm wax-like consistency. It contains acrylates, waxy components, fillers as well as light activators that keep polymerisation contraction extremely low.

### Specific precautions

When modelling, work by heating under a fume cupboard or with good ventilation.

### Application

**How is LiWa cured?** LiWa is cured by light polymerisation. The light activators are set to activate curing from UV light to the visible range. So in the range between 280 - 520 nm. Any standard laboratory light-curing device can be used for this purpose, whether with halogen or UV lamp or with a combination of both. It is important to ensure that the object is sufficiently cooled during polymerisation.

**How long does it take for LiWa to cure?** The curing time depends on the type of lamp, power of the lamp, age of the lamp and the colour of the plaster model. A dark plaster colour reflects less light than a light plaster colour. Fluorescent components in the plaster shorten the polymerisation time. For example, a LiWa thickness of 2-3 mm in the Uni XS® device from Heraeus-Kulzer takes an average of 15 - 20 seconds. (From the fastest light unit e.g. Uni XS 30 sec. to normal unit approx. 2 - 4min).

**How is the correct polymerisation time determined?** The following test has proved successful.

Apply a strand of LiWa material about 3 cm long and 2-3 mm in diameter to an isolated laboratory plaster model. The bottom part of the strand should be well adapted to the insulated model. After a pre-determined polymerisation time (e.g. 3 minutes), the curing is checked by detaching the strand from the model and cutting off approx. 2-3 mm with a thin cut-off wheel. The separated part must be thoroughly cured. If there are still soft areas in the inner core or on the model side, the polymerisation time must be extended and a cut test must be carried out again until complete curing has taken place! Please note! Light type, intensity, time, plaster colour, LiWa layer thickness and the age of the lamp play an important role. Therefore, regularly carry out test polymerisations as described above as part of quality assurance.

**How does LiWa behave after curing?** After light polymerisation, LiWa can be processed mechanically like plastic with rotating tools such as milling cutters, diamond grinders, sandpaper and rubber polishers. LiWa is stable enough not to deform and at the same time flexible enough, for example, to lift a clamp easily even from underneath itself. The resilience is excellent and a particular strength of the new material.

**How is LiWa processed?** LiWa can be modelled "hot" and "cold" (by kneading and rolling). For hot modelling, the easiest way to melt is with an electric wax knife at a temperature of 150 to 250°C. If the absorbed material already "steams" when the wax knife is dipped into the mass, the temperature should be reduced step by step. If it sticks to the electric wax knife, the temperature must be increased. LiWa behaves like a lead-dead wax, i.e. it is non-resetting and does not contract. When melting, LiWa flows into a thin layer, resulting in a homogeneous smooth surface that does not change perceptibly during and after polymerisation. Since LiWa remains viscous for longer when working hot, it can be cooled down to a slightly better workable consistency with a cold spray (LiWa-Cool). In cold modelling (i.e. at room temperature), LiWa can be plastically moulded. The consistency is roughly comparable to modelling clay. The desired plastic consistency is achieved by rolling with the fingers. The prerequisite is the wearing of protective gloves, as prescribed by the Employer's Liability Insurance Association when working with acrylates and acrylic-containing materials, in order to exclude a possible sensitisation reaction or an allergic reaction. In order to prevent the modelled sections from crumpling, fingers should be used with water. As a further advantage, this also results in cooling of the modelling material. As is usual in the resin veneering technique, metal and silicone tip instruments are also very helpful in cold modelling. In addition, the basic mould should be slightly oversized for cold modelling and shaped with rotating instruments after light polymerisation. Small additions or the finishing of the fine modelling can be done quickly and normally with the usual modelling waxes.

**How is plaster or plastic insulated against LiWa?** Plaster insulation: With LiWa Iso Step I, the plaster surface is sealed without applying (diffuses into the plaster pores). Application time min. 50 seconds. Then insulate with LiWa Iso Step II (blow away excess). Application time min. 30 seconds. Plastic insulation: With the new LiWa universal insulator, the model is insulated 2x in a row. Application time 30 seconds each. Can also be used for insulation LiWa against LiWa, e.g. with attachment parts, wax, plastics and metals. Important! Model sections/ dies must not be blocked out with a light-curing block-out material and likewise not be varnished with a light-curing die varnish. Despite careful isolation, an inseparable bond with the light-curing block-out material or the die spacer may occur after polymerisation of LiWa.

**How is LiWa blocked out?** The necessary blocking out for model casting can be done as usual with tin foil or wax, which can be adapted exactly to the course of the jaw. It must also be pre-treated with LiWa Universal Isolant 1x.

**How can LiWa parts be connected to each other?** Uncured LiWa parts can be joined easily by smearing them with the electric wax knife. Already polymerised and ground parts must be activated with LiWa Akti Bond at the bonding point.

**How can the inhibition layer (light lubricating film) be removed?** After light polymerisation, LiWa forms an inhibition layer on the surface in combination with air/oxygen like any other light material or light material combination. This lubricating film can be easily removed with alcohol or monomer.

**How can the finished LiWa surface be tempered/smoothed?** By applying a thin layer of LiWa Finish-Varnish (this is a crystal-clear surface sealing lacquer), you obtain a mirror-smooth compact surface. On the one hand, this facilitates desvesting and, on the other hand, it considerably shortens the working time in the metal. Likewise, the compact LiWa Finish Varnish layer gives additional stability to a very graceful work made of LiWa. LiWa Finish Varnish is also very suitable for sealing the surface of wax models, for reinforcing contact points and as a retention adhesive. Important! In the case of large areas/parts to be sealed, however, LiWa Finish Varnish may only be applied partially, as otherwise (in the case of large surface treatment) the light polymerisation may cause the entire work to distort. Due to the transparency of the varnish, the polymerisation time is kept to a minimum.

**How is LiWa applied?** LiWa can be used with the 2 or 4 mm round wires that are part of the system. However, conventional wax wires can be used just as well. Here, the known and generally valid rules for the attachment point of the wires and for the casting of the wax objects (depending on the casting equipment) must be observed.

**How is LiWa embedded?** LiWa is basically a plastic material. It burns without residue during preheating. Like all plastics, LiWa also tends to swell during the heating phase. Therefore, the choice of investment material is of particular importance. It should be stable and not tear even when LiWa swells. With controlled heating, these problems occur very rarely. With speed casting, however, you should first make a few test castings before you have found the optimum mass. By dispensing with an investment model, not only model casting investments but also conventional K&B investments are suitable. These have the advantage of an even finer surface. Very good experiences have been made with Heravest Speed (Heraeus-Kulzer), Formula fast (Dentona AG) and Vulcan (HM Dentalbedarf), for example.

**Disclaimer of warranty:** The above information has been tested in numerous laboratory and practical trials. The user is responsible for correct application and processing. Errors occurring at the user's premises as a result of the work do not constitute grounds for a warranty claim.

### Warnings

#### Residual risks

The product has only an indirect clinical use, as it is applied to make metal parts (crowns, bridges) in the laboratory.

#### Side effects

See "Contraindications" and "Special precautions".

#### Interactions

No interaction of the product and its metabolites with other medical devices, medicinal products or other substances is known of.

#### Storage

At elevated room temperatures (hot summer days >27°), LiWa should be stored in a cool and dark place. This makes the processing time longer and prevents the material from "sticking". Always keep LiWa in the original opaque box/packaging. Only as much LiWa material should be removed as is actually needed. LiWa material must not be exposed to natural or artificial light for a long time, otherwise it will polymerise.

#### Minimum shelf life

Do not use after the expiry date.

#### Variants

<a href="#">REF</a> WP5075	LiWa, 1 x 20 g
<a href="#">REF</a> WP5073	LiWa II, 1 x 20 g
<a href="#">REF</a> WP5090	LiWa Form grained plates, 12 pcs.
<a href="#">REF</a> WP5097	LiWa Form smooth plates, 12 pcs.
<a href="#">REF</a> WP5093	LiWa Form molar clamps, 10 x 18 pcs.
<a href="#">REF</a> WP5094	LiWa Form round hole retentions, 12 pcs.
<a href="#">REF</a> WP5092	LiWa Form UK bracket, medium, approx. 4mm, 30 pcs.
<a href="#">REF</a> WP5096	LiWa Form wax wire, 2 mm, 30 pcs.
<a href="#">REF</a> WP5096A	LiWa Form wax wire, 4 mm, 20 pcs.
<a href="#">REF</a> WP5096B	LiWa Form wax wire, 1 mm, 30 pcs.
<a href="#">REF</a> WP5076	LiWa Form retention grid, round holes, 10 pcs.
<a href="#">REF</a> WP5089	LiWa Form retention grid, square holes, 10 pcs.
<a href="#">REF</a> WP6135	LiWa Form grained plates, large size, approx. 150x75mm, 10 pcs.
<a href="#">REF</a> WP6129	LiWa Dip dipping wax, transparent, 2 x 40g in tubes
<a href="#">REF</a> WP5140	LiWa Flow modelling resin, 1 x 2 ml in syringe
<a href="#">REF</a> WP5079	LiWa COOL cold spray, 1 x 200 ml
<a href="#">REF</a> WP5095	LiWa ISO Step I + II, 2 x 10 ml
<a href="#">REF</a> WP5082	LiWa Akti Bond, 1 x 10 ml
<a href="#">REF</a> WP5074	LiWa Finish Varnish, 1 x 10 ml
<a href="#">REF</a> WP5083	LiWa Universal-Isolant, 1 x 10 ml
<a href="#">REF</a> WP5100X	LiWaxer Set
<a href="#">REF</a> WP5100NX	LiWaxer NX Set
<a href="#">REF</a> WP6133	LiWa Dip Set

#### Disposal

Uncured material can be cured and disposed of with household waste. Uncured material should be disposed of like methacrylate waste.

**All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State where the user and/or the patient is established.**

We know how

**dent a pharm**  
Produktionsgesellschaft mbH  
Schusterring 35  
25355 Barmstedt/Germany  
Fon: +(49) 4123/9225-0  
Fax: +(49) 4123/9225-48  
www.dent-a-pharm.de  
info@dent-a-pharm.de

**Willmann & Pein GmbH**  
Schusterring 35  
25355 Barmstedt/Germany  
Fon: +(49) 4123/9228-0  
Fax: +(49) 4123/9228-49  
www.wp-dental.de  
info@wp-dental.de

H315 causes skin irritation  
H317 may cause allergic skin reactions  
H319 causes serious eye irritation  
H400 very toxic to aquatic organisms

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