# Aquafluor®

- Seite 2 - de: Gebrauchsanweisung  
Page 3 - en: Instructions for use  
Page 4 - fr: Manuel d'instruction  
Página 5 - es: Instrucciones de uso  
Pagina 6 - it: Istruzioni per l'uso  
Strana 7 - bs: Uputstvo za upotrebu  
Pagină 8 - ro: Instrucțiuni de utilizare  
Sayfa 9 - tr: Kullanma kılavuzu  
Strona 10 - pl: Instrukcja użytkowania  
Стр. 11 - ru: Инструкция по применению

*We know how*

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**Patientenzielgruppe**

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Zahnarzt

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Qualitative Angaben:

Natriumfluorid, Bis-(hydroxyethyl)-aminpropyl-N-hydroxyethyl-octadecylamin-dihydrofluorid, wässrige Methacrylat-Dispersion

Quantitative Angaben:

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**Composition**

Qualitative data:

Sodium fluoride, bis-(hydroxyethyl)-aminpropyl-N-hydroxyethyl-octadecylamine dihydrofluoride, aqueous methacrylate dispersion

Quantitative data:

Sodium fluoride, bis-(hydroxyethyl)-aminpropyl-N-hydroxyethyl-octadecylamine dihydrofluoride 0.15 % calculated as pure fluoride

**Intended purpose of the product**

Fluoridation of the tooth enamel.

**Indications**

Aquafluor is a varnish for fluoridation of the tooth enamel and desensitisation of sensitive teeth.

**Contraindications**

Do not use in patients with a known allergy to methacrylates.

**Performance characteristics of the product**

After varnish formation (application of dispersion and dry blowing) on the enamel, the varnish softens and releases fluoride.

**Specific precautions**

None

**Properties**

Aquafluor is an air-drying dispersion that forms a well-adhering film on enamel and dentin. The adherent film softens when it comes into contact with saliva or water for a longer time. In the process, the fluorides contained in the polymer matrix can diffuse to the dentin or enamel surface and exert their effect there. Due to this behaviour, minimal amounts of fluoride are completely sufficient for optimal fluoridation. Once the varnish is completely soaked, the fluoride release is complete. The varnish matrix then detaches from the enamel or dentin surface.

**Application**

Clean surfaces to be treated and, if possible, dry them. Transfer Aquafluor drop by drop to the brush (it is not necessary to shake up the aqueous dispersion before removal) and touch surfaces thinly. Dry the dispersion carefully and extensively, preferably with warm air, until the polymer film has completely formed and is dry. Do not brush teeth for about 24 hours after treatment. When used correctly, the polymer layer will not come off until it is completely soaked and maximum fluoride release has occurred.

**NOTE**

The application aids included with the product (cannulas, disposable brushes, etc.) are only for single use on a single patient!

**Duration of application**

For caries prophylaxis, treatment should be carried out 2–3 times a year (or 3–4 times if there is an increased susceptibility to caries).

**Warnings****Residual risks**

The processing instructions and precautionary measures described in the instructions for use of the product must be observed in every instance. If, despite following these instructions, an adverse reaction occurs in the patient or other adverse events occur, please inform us in every case, describing the accompanying circumstances as precisely as possible.

**Side effects**

Undesirable side effects of this medical product are extremely rare, provided the product is processed and used properly. If you become aware of any adverse effects – even in doubtful cases – please inform us and describe the accompanying circumstances and symptoms as precisely as possible. We follow up on every indication.

No risks in the event of an overdose are known of.

**Interactions**

Interaction of the product with other medical devices, medicinal products or other substances is not known.

**Storage**

The Product should be stored between 0° C - 22°C.

**Minimum shelf life**

Do not use after the expiry date.

**Variants**

**REF**WP4005 1 x bottle with 5.0 ml liquid + 50 brushes

**REF**WP4004 1 x bottle with 10.0 ml liquid

**REF**WP4949 1 x bottle with 5.0 ml liquid, Strawberry + 50 brushes

**REF**WP4948 1 x bottle with 10.0 ml liquid, Strawberry

**Disposal**

Dilute unused liquid with water and dispose of with household waste. Cured residues can also be disposed of with household waste.

**All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient is established.**

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**Willmann & Pein GmbH**