# C-Bond Self-Etch®

|           |       |                           |
|-----------|-------|---------------------------|
| Seite 2   | - de: | Gebrauchsanweisung        |
| Page 3    | - en: | Instructions for use      |
| Page 4    | - fr: | Manuel d'instruction      |
| Página 5  | - es: | Instrucciones de uso      |
| Pagina 6  | - it: | Istruzioni per l'uso      |
| Strana 7  | - bs: | Uputstvo za upotrebu      |
| Pagină 8  | - ro: | Instrucțiuni de utilizare |
| Sayfa 9   | - tr: | Kullanma kılavuzu         |
| Strona 10 | - pl: | Instrukcja użytkowania    |
| Стр. 11   | - ru: | Инструкция по применению  |

*We know how*

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**Patientenzielgruppe**

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## Instructions for use

**Patient target group**

Children, teenagers, adults and senior citizens

**Intended user**

Dentist

**Composition**

Qualitative data:

Oligomeric methacrylates, HEMA, camphorquinone

Quantitative data in % G/G:

Methacrylates 99.5 Additives 0.5

**Intended purpose of the product**

C-Bond Self-Etch is a self-etching, light-curing, one-component bonding agent for the adhesive restoration technique.

**Indications**

Bonding agent for the adhesive restoration technique of dental hard tissue defects

- Bonding of composite materials and compomers to enamel and dentine
- Adhesive bonding of composite restorations on enamel and dentine

**Contraindications**

Product contains methacrylates and organic acids. The product should not be used in the case of known hypersensitivities and allergies to these ingredients.

Do not apply directly or indirectly to the pulp (perform direct or indirect pulp capping beforehand).

**Performance characteristics of the product**

Adhesive bonding of resin fillings to enamel and dentine

**Specific precautions**

None

**Application**

1. Conventional filling technology

## a) Cavity preparation

According to the principle of adhesive restoration technique.

If necessary, a punctiform capping of the pulp can be performed with a hard-setting calcium hydroxide cement. Processing of enamel with a diamond grinder. For unpolished enamel, etch the entire cavity with etching gel starting from the enamel margin (30 s for enamel, max. 15 s for dentin). Then siphon off the etching gel, rinse for approx. 20 s and remove excess moisture with an air blower. Do not overdry the dentin. The aim is to achieve a slightly moist dentin surface. If the etched material does not appear matt/chalky white, repeat the etching process. The use of a rubber dam is recommended for hygienic reasons.

## b) Application

Shake the bottle well before use and add a drop to the dot mould. Massage well with a disposable brush C-Bond Self-Etch 20 s on dentin and enamel. Gently blow the C-Bond Self-Etch 5 s with oil-free compressed air. Then light-cure for at least 20 s (halogen or LED polymerisation units with an illumination intensity  $\geq 600 \text{ mW/cm}^2$ ). With lower light output, light curing is prolonged accordingly.

## c) Composite and compomer fillings

Apply restoration material according to the respective manufacturer's instructions.

**Processing instructions:**

The working instructions and precautionary measures described in the package leaflet accompanying the product must be observed in every instance. Should an undesirable reaction occur in the patient despite compliance or other undesirable events occur, we ask in any case for notification with as precise a description of the accompanying circumstances as possible.

**Precautions**

The user and patient should wear protective goggles before application, as well as when applying etching gel.

Avoid eye contact.

In case of eye contact, rinse immediately with plenty of water and seek medical advice.

Avoid contact with the skin.

In the event of skin contact, remove with a cloth and alcohol and then wash thoroughly with soap and water. In case of hypersensitivity of the skin or skin rash, discontinue use and consult a doctor.

Details

Delivery only to dentists and dental technicians or on their behalf. Keep out of the reach of children! Only completely emptied packaging should be recycled.

**Side effects**

Undesirable side effects of this medical device are not to be expected if used correctly. If you become aware of any adverse effects – even in doubtful cases – please inform us and describe the accompanying circumstances and symptoms as precisely as possible.

**Interactions**

Phenol-containing substances or other substances (e.g. zinc oxide eugenol) that could influence the polymerisation must not be brought into contact with the C-Bond Self-Etch.

**Storage**

Do not expose C-Bond Self Etch to direct sunlight or other heat sources.

**Minimum shelf life**

3 years

Do not use after the expiry date.

**Variants**

**REF**WP4048 C-Bond Self-Etch, 5.0 ml

**REF**WP4018 C-Bond Self-Etch, 10.0 ml

**Disposal**

Unused material can be cured and disposed of with household waste. Uncured material should be disposed of like methacrylate waste.

**Any serious incident that has occurred in relation to the product must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.**

We know how

**dent a pharm**  
Produktionsgesellschaft mbH  
Schusterring 35  
25355 Barmstedt/Germany  
Fon: +(49) 4123/9225-0  
Fax: +(49) 4123/9225-48  
www.dent-a-pharm.de  
info@dent-a-pharm.de

**Willmann & Pein GmbH**  
Schusterring 35  
25355 Barmstedt/Germany  
Fon: +(49) 4123/9228-0  
Fax: +(49) 4123/9228-49  
www.wp-dental.de  
info@wp-dental.de

Made in Germany

H317 may cause allergic skin reactions

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