Pro-Fill®

- Seite 2 - de: Gebrauchsanweisung
Page 3 - en: Instructions for use
Page 4 - fr: Manuel d'instruction
Página 5 - es: Instrucciones de uso
Pagina 6 - it: Istruzioni per l'uso
Strana 7 - bs: Uputstvo za upotrebu
Pagină 8 - ro: Instrucțiuni de utilizare
Sayfa 9 - tr: Kullanma kılavuzu
Стр. 10 - ru: Инструкция по применению

We know how

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Patientenzielgruppe

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Instructions for use

Patient target group

Children, adolescents, adults and older patients

Intended users

Dentist

Composition

Qualitative and quantitative data:

1 g contains 229 mg of silicon dioxide, 595 mg of urethane methacrylate, 176 mg of TEDMA

Intended purpose of the product

Temporary closure of cavities.

Indications

Suitable for temporary fillings, temporary protection involving inlay and onlay techniques.

Contraindications

Do not use in patients with known allergy to acrylates.

Performance characteristics of the product

Pro-Fill is a temporary closure material that cures with light.

Specific precautions

None

Application

1. Prepare, underfill and take the impression as usual, rinse thoroughly with water and dry the cavity well.
2. Pro-Fill is placed in the cavity with a suitable instrument (flat plastic spatula, spherical stopper, etc.).
3. Smooth the surface with the help of a matrix.
4. The temporary filling is then cured with a halogen lamp for at least 20-40 seconds. The curing time depends on the layer thickness of the filling (e.g. 1 mm = 5 sec.; 5 mm = 20 sec.)
If the halogen lamp cannot be brought close enough to the place to be cured, the curing time must be extended. Light devices should emit or emit light with a wavelength of 450 nm. The intensity of radiation should be at 500 mW per cm² be.
5. Check the surface. Excess material can be removed with a scalpel. Work and finish the provisional filling as usual.

Removing the temporary filling

Pressing a pointed instrument (probe) into the temporary filling allows it to be removed easily. The preparation is not damaged, as Pro-Fill remains elastic in the cured state.

Warnings

Residual risks

The working instructions and precautionary measures described in the instructions for use of the product must be observed in all cases. If, despite following these instructions, an adverse reaction occurs in the patient or other adverse events occur, we request that you inform us in every case, describing the accompanying circumstances as precisely as possible.

Side effects

Undesirable side effects of this medical device, when the device is properly processed and used, are extremely rare. If you are aware of undesirable side effects – even in doubtful cases – we ask in any case for notification with as precise a description of the accompanying circumstances and symptoms as possible. We follow up on every indication.

No risks in the event of an overdose are known of.

Interactions

Interaction of the product with other medical devices, medicinal products or other substances is not known.

Storage

The Product should be stored between 0° C - 22°C. Storage in the refrigerator is possible. Bring to room temperature before use. Store away from light!

Minimum shelf life

Do not use after the expiry date.

Variants

 WP4041 2 x 4.5 g in screw syringes

Disposal

Unused material can be cured and disposed of with household waste. Uncured material should be disposed of like methacrylate waste.

All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient is established.

We know how

dent a pharm
Produktionsgesellschaft mbH
Schusterring 35
25355 Barmstedt/Germany
Fon: +(49) 4123/9225-0
Fax: +(49) 4123/9225-48
www.dent-a-pharm.de
info@dent-a-pharm.de

Willmann & Pein GmbH
Schusterring 35
25355 Barmstedt/Germany
Fon: +(49) 4123/9228-0
Fax: +(49) 4123/9228-49
www.wp-dental.de
info@wp-dental.de

Made in Germany

H317 may cause allergic skin reactions

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