Securafil®

Seite 2	- de: Gebrauchsanweisung
Page 3	- en: Instructions for use
Page 4	- fr: Manuel d'instruction
Página 5	- es: Instrucciones de uso
Pagina 6	- it: Istruzioni per l'uso
Strana 7	- bs: Uputstvo za upotrebu
Pagină 8	- ro: Instrucțiuni de utilizare
Sayfa 9	- tr: Kullanma kılavuzu
Strona 10	- pl: Instrukcja użytkowania
Стр. 11	- ru: Инструкция по применению

We know how

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Patientenzielgruppe

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Instructions for use

Patient target group

Children, adolescents, adults and older patients

Intended user

Dentist

Composition

Qualitative composition:

Powder: Calcium aluminium fluorosilicate, polyacrylic acid, iron oxide

Liquid: Polyacrylic acid, tartaric acid, water

Quantitative data:

Powder: > 90 % calcium-aluminium-fluorosilicate

Liquid: > 40 % polyacrylic acid, > 50 % water

Intended purpose of the product

Securafil is a glass-ionomer cement preferred for filling Class III and V cavities.

Indications

For filling class III and V cavities.

Contraindications

Do not use in case of allergies to the ingredients.

Performance characteristics of the product

The cement hardens and chemically adheres to the dentin and enamel.

Adhesion to dentin and enamel occurs via the reaction of carboxylate groups of the polycarboxylic acid with the calcium of the hydroxy-apatite and the collagen. This direct chemical bond allows much more healthy tooth material to be preserved, as no undercuts are necessary.

Specific precautions

Do not allow the unbound material, the tooth cleaner or the varnish to come into contact with the eyes. If such contact should occur accidentally, wash the eyes immediately with plenty of water.

Application

Preparation of the cavities:

Healthy tooth material should be preserved as much as possible, as no undercuts are required for better adhesion. The edge should not be thinner than 0.5 mm.

Pulp protection:

Rinse the prepared cavities well with water and dry them in an oil-free air stream. For very deep cavities, use calcium hydroxide preparations (e.g. Calcident 450) as pulp capping.

Preparation of the dentine:

Apply a tooth cleaner to the freshly prepared dentin surfaces and to remove the plaque (allow to act for 10 seconds here), rinse with water and air dry. Avoid dehydration!

Mixing Securafil:

Shake the powder well before removing it! Two level measuring spoons are used for two full drops of liquid (holding the bottle vertically without air bubbles, squeeze out drops with slight pressure). Mixing is done first with half of the powder, then further with half of the remaining powder and so on until the desired consistency is reached. To minimise evaporation of the water, only a small area of the mixing block should be used.

Application of Securafil:

The mixed material is placed directly into the cavity, where it is tamped down and, if necessary, the surface is shaped with a suitable matrix or coated with a varnish (e.g. Secura Varnish). Securafil is workable for 2 minutes (including mixing time). Once the curing process begins, the material should no longer be manipulated. Curing time is approx. 3 minutes.

Trimming:

After curing, the matrix can be removed. The now unprotected part should be covered with varnish. After another 6 minutes, the trimming should begin. The excess material is removed using precious corundum grinders or flexible grinding wheels coated with Vaseline. A complete polishing is done with Vaseline-coated grinding stones or discs. Attention! Do not use water spray! Afterwards, the freshly prepared area is again covered with a varnish to protect the ionomer cement from water during the first 24 hours.

Warnings Residual risks

The working instructions and precautionary measures described in the package leaflet accompanying the product must be observed in any case. If an adverse reaction occurs in the patient or other adverse events occur (despite following these instructions), we request that you inform us in every case, describing the accompanying circumstances as precisely as possible.

Note: Clean instruments under water immediately after use.

Side effects

Undesirable side effects of this medical device, when the device is properly processed and used, are extremely rare. If you become aware of any adverse effects – even in doubtful cases – please inform us and describe the accompanying circumstances and symptoms as precisely as possible.

Interactions

Interaction of the product with other medical devices, medicinal products or other substances is not known.

Storage

Store in a dry place. The Product should be stored between 0° C - 22°C.

Minimum shelf life

Securafil has a shelf life of 3 years if unopened and 6 months if opened.

Variants

 WP4142 - Securafil, A2, 15.0 g powder / 10.0 g liquid + measuring spoon

 WP4143 - Securafil, A3, 15.0 g powder / 10.0 g liquid + measuring spoon

 WP4095 - Securafil, A2, 7.0 g powder / 5.0 g liquid + measuring spoon

 WP4096 - Securafil, A3, 7.0 g powder / 5.0 g liquid + measuring spoon

 WP4145 - Securafil, 10.0 g liquid

Disposal

Dispose of unused material according to regulations (hazardous waste).

All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State where the user and/or the patient is established.

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Willmann & Pein GmbH